

**Oxxiús**  
Simply Light

# **LaserBoxx**

## **Single Frequency Lasers**



Raman Spectroscopy  
Brillouin Scattering  
Interferometry  
Dynamic Light Scattering  
Holography  
Laser Doppler Velocimetry  
Shearography

# LaserBoxx


One platform for all colors

The single frequency lasers from Oxixus have been designed with versatility in mind. They utilize advanced technologies such as DPSS lasers and stabilized laser diode.

## Technology

### DPSS lasers


The LCX and LPX *LaserBoxx* are diode-pumped solid-state (DPSS) laser sources. The unique feature of these models is a proprietary, **Alignement-free Monolithic Resonator (AMR)**.



The elements of the resonator are assembled into a single ultra-low-loss optical subsystem, using a proprietary crystal bonding technique.

#### Benefits of the AMR

This technology yields to highly efficient pump schemes and allows for the highest spectral quality and wavelength stability ( $\leq 1\text{pm}$ ) on the market, as well as an important robustness over time. The LCX and LPX models are insensitive to temperature variations and to mechanical vibrations.



### Diode lasers

The LBX and LSX models are based on integrating a **temperature-stabilized laser diode**.

#### Benefits of VBG stabilized lasers

These models deliver an ultra-narrow linewidth emission due to their stable design and their proprietary wavelength-locking routine.

#### Common key features

- Single frequency
- Narrow linewidth
- Up to 500 mW continuous wave
- Integrated control electronics
- Low profile laser head
- SM/PM/MM fiber coupling options
- USB and RS-232 interfaces
- 100 x 40 mm<sup>2</sup> Industry standard footprint (LBX and LCX)


#### Common specifications

- Power stability (over 8 h and  $\pm 3\text{ K}$ )  $\pm 1\%$
- Power adjustment optional with L1C-MPA/AOM
- Optical noise (10 Hz - 20 MHz bandwidth)  $\leq 0.2\%$ , 2 m

## Performances


Beam profile

LBX-633



Narrow linewidth spectrum

LSX-785S-ISO



# Optical specifications

## DPSS lasers


	Emission wavelength (Tolerance)	Optical power (Continuous wave)	Wavelength stability over 8 h and $\pm 3$ K	Linewidth	Side mode suppression ratio	Polarization state	Beam waist diameter (at $1/e^2$ )	Beam quality factor ( $M^2$ )	Beam circularity in far field	Coupling efficiency Into a single mode fiber				
LCX-532S	532.3 nm ( $\pm 0.5$ nm)	50 / 100 / 150 200 / 300 mW	$\leq 1$ pm	$\leq 1$ MHz	$\geq 30$ dB	linear, vertical 100:1 extinction ratio	0.7 mm ( $\pm 0.1$ mm)	$\leq 1.1$	$\geq 90\%$	$\geq 70\%$				
LPX-532S		500 / 800 mW												
LCX-553S	553.0 nm ( $\pm 0.4$ nm)	50 / 100 200 mW												
LCX-561S	561.4 nm ( $\pm 0.4$ nm)	50 / 100 150 / 200 mW												
LPX-561S		300 / 500 mW												
LCX-946S	946.0 nm ( $\pm 0.3$ nm)	50 mW				linear, vertical 300:1 extinction ratio								
LCX-1064S	1064.6 nm <sup>(1)</sup> ( $\pm 0.6$ nm)	100 / 200 / 300 500 mW												

<sup>(1)</sup> The LCX-1064S also emits a class 1 visible aiming beam.


## Diode lasers

	Emission wavelength (Tolerance)	Optical power (Continuous wave)	Wavelength stability over 8 h and $\pm 3$ K	Linewidth	Side mode suppression ratio	Linear polarization extinction ratio	Beam waist diameter (at $1/e^2$ )	Beam quality factor ( $M^2$ )	Beam circularity in far field	Coupling efficiency into a single mode fiber
LBX-633S	633 nm ( $\pm 0.5$ nm)	40 mW	$\leq 10$ pm	$\leq 100$ MHz	$\geq 35$ dB	50:1	0.2 to 0.6 mm	$\leq 1.9$	$\geq 65\%$	$\geq 50\%$
LSX-785S-ISO	785 nm ( $\pm 0.5$ nm)	150 mW with isolator			$\geq 25$ dB	100:1	0.4 to 0.6 mm	$\leq 1.25$	$\geq 90\%$	
LBX-830S	830 nm ( $\pm 0.5$ nm)	150 mW			$\geq 35$ dB	50:1	0.5 to 1.0 mm	$\leq 1.9$	$\geq 65\%$	
LBX-785S-MM	785 nm ( $\pm 0.5$ nm)	450 / 600 mW	$\leq 10$ pm	100 pm 60 pm typ.	$\geq 35$ dB	Delivery on a multimode fiber 100 $\mu$ m-core diameter 0.22 numerical aperture				
LBX-830S-MM	830 nm ( $\pm 0.5$ nm)									
LBX-976S-MM	976 nm ( $\pm 0.5$ nm)									
LBX-1064S-MM	1064 nm ( $\pm 0.5$ nm)									


Power stability versus temperature  
LCX-561S



Wavelength stability versus time  
LCX-561S



Narrow linewidth spectrum  
LCX-532S





## L1C

The L1C platform offers an efficient, compact and cost-effective solution to add advanced features to the **LaserBoxx** lasers:

### Motorized Power Attenuator (MPA)

Power adjustment; the MPA is a continuous control of the optical power. It is ideal to ensure wavelength stability.

- 0-100 % power range (full contrast with e-m shutter)
- Analog voltage or software command inputs
- Transmission ratio:  $\geq 85$  %
- Response time:  $< 1$  s
- Dynamic range  $\geq 30$  dB
- Digital modulation up to 5 Hz with e-m shutter
- Compatible with an isolator

### Acousto-Optic Modulator (AOM)

Modulated output; analog or digital modulation up to 3 MHz.

- Transmission ratio:  $\geq 85$  %
- Response time: 100 ns typical
- Dynamic range  $\geq 40$  dB
- Electro-mechanical shutter for a complete extinction
- Plug & play version with BTC-AOM driver
- USB and Ethernet interfaces

### Isolator (ISO)

An optical isolator protects the laser source in reflective environments.


- Degree of isolation  $\geq 20$  dB typ.
- Transmission ratio  $\geq 70$  % typ. (wavelength-dependent)
- Compatible with the MPA



**L1C-532S**

## Fiber coupling

A rugged and compact accessory that injects the laser beam into a single mode (SM) fiber, a polarization maintaining (PM) fiber, or a multimode (MM) fiber.



**Fiber coupling**


	SM and PM Fiber	MM Fiber (50 $\mu$ m, 0.22 NA)
Coupling efficiency	LCX LPX LSX $\geq 70$ % LBX-S $\geq 50$ %	$\geq 80$ %
Power stability over 8 h, $\pm 1.5$ K	$\pm 2$ %	$\pm 2$ %
Polarization extinction ratio (PMF only)	$\geq 50:1$	n/a
Available optical connectors	FC-APC	AR-coated SMA FC-APC
Fiber length	2.0 m	2.0 m

# System specifications


## Plug and Play, CDRH-compliant versions

Dimensions in mm

**LCX series**  
DPSS laser




**LPX series**  
DPSS laser




The electro-mechanical shutter is standard on LPX models

**LBX series**  
Stabilized laser diode



Optional heatsink available


**LSX series**  
Stabilized laser diode



Optional heatsink available


### Power-adjustable versions (PPA)

**PPA - ControlBoxx**




### Fixed power versions (PPF)

**PPF - RemoteBoxx**



### For improved stability

**Heatsink**




<b>Compliance</b>	CE and FDA 21 CFR 1040.10 / 1040.11
<b>Operating temperature</b>	10 to 38 °C (ambient air)
<b>Power consumption</b>	≤ 25 W
<b>Storage temperature</b>	0 to 60 °C

<b>Supply voltage</b>	100 to 240 VAC external power supply
<b>Warm-up time</b>	LCX, LPX: ≤ 10 minutes LBX, LSX: ≤ 2 minutes
<b>Communication interfaces</b>	USB, RS-232, dedicated interface


## OEM - dedicated versions

Control electronics is integrated into the laser head

**LCX series**  
DPSS Laser




**LPX series**  
DPSS Laser




The electro-mechanical shutter is standard on LPX

**LBX series**  
Stabilized laser diode



**LSX series**  
Stabilized laser diode  
Integrated isolator



	LCX and LPX	LBX and LSX
<b>Compliance</b>	FDA 21 CFR 1040.10 / 1040.11	
<b>Operating temperature baseplate</b>	10 to 50 °C	20 to 35 °C
<b>Power consumption</b>	≤ 25 W	≤ 10 W
<b>Storage temperature</b>	0 to 60 °C	
<b>Supply voltage</b>	5 to 12 VDC	
<b>Warm-up time</b>	≤ 10 minutes	≤ 2 minutes
<b>Communication interfaces</b>	USB, RS-232, dedicated electronical interface	



## Our sales network

Our distributors are present all over the world, making our products easily accessible wherever you are. To find the full list of our partners and their locations, visit our website: [www.oxxius.com/contact-us](http://www.oxxius.com/contact-us).



## About Oxxius

Founded in 2002, Oxxius develops, manufactures and sells advanced DPSS and laser diode modules across the ultraviolet, visible, and near-infrared spectra.

Our solutions deliver exceptional optical performance in an ultra-compact design, making them easy to integrate into instruments for life science, metrology, and manufacturing applications.

Oxxius also offers compact and versatile multicolor laser sources wavelength combiners, with up to 7 laser lines.



OXXIUS  
4 rue Louis de Broglie - 22300 Lannion, France  
Phone: +33 296 48 70 28  
E-mail: [sales@oxxius.com](mailto:sales@oxxius.com)  
[www.oxxius.com](http://www.oxxius.com)

