

# LCX-532S

DPSS laser

## Optical characteristics \*

Emission wavelength	532.3 nm ( $\pm 0.3$ nm tolerance)
Spectral linewidth, FWHM	$\leq 1$ MHz
Wavelength stability over 8 hours, temperature within $\pm 3^\circ\text{C}$	$\leq 1$ pm
Coherence length	$\geq 100$ m
Side mode suppression ratio	$\geq 30$ dB

## — Optical power

Part number	Free space output	Fiber-coupled output	
		Single mode	Multi mode fiber
LBX-532S-50-CSB-(OE/PPF)	50 mW	35mW	40 mW
LBX-532S-100-CSB-(OE/PPF)	100 mW	70 mW	80 mW
LBX-532S-150-CSB-(OE/PPF)	150 mW	105 mW	120 mW
LBX-532S-200-CSB-(OE/PPF)	200 mW	140 mW	160 mW

**Power stability**  $\pm 1\%$   
over 8 hour, temperature within  $\pm 3^\circ\text{C}$

**Control modes** Automatic Current Control (ACC)

**Optical noise**  $\leq 0.2\%$  rms  
10 Hz-20 MHz bandwidth

## — External modulation

Analog Modulation	L1C-MPA option	L1C-AOM option
Rise / fall time	$\leq 1$ s	$\leq 300$ ns
-3dB modulation bandwidth	-	$\geq 3$ MHz
Digital Modulation		
Rise / fall time	$\leq 30$ ms	$\leq 300$ ns
-3dB modulation bandwidth	-	$\geq 3$ MHz

## — Transverse singlemode free-space beam

**Beam diameter (typ)** at  $1/e^2$ ,  
5 cm from the beam aperture 0.7 mm

**Beam divergence** at  $1/e^2$   
full beam, in far field 1 mrad

**Beam quality factor**  $M^2 \leq 1.1$


**Beam circularity**, in far field  $\geq 90\%$

**Polarization state** Linear, vertical at  $\pm 5^\circ$


**Polarization**  $> 1000 : 1$

## General specifications

	Plug and Play version	OEM version
Compliance	CE FDA 21 CFR 1040.10/1040.1	FDA 21 CFR 1040.10 / 1040.11
Operating temperature	10 to 38°C ambient air with optional heat sink	20 to 35°C baseplate
Power consumption	20 W max.	
Storage temperature	0 to 60°C	
Humidity	Up to 90%RH, non-condensing	
Supply voltage	100 to 240 VAC external power supply	5 to 12 VDC
Warm-up time	$\leq 2$ minutes	
Interfaces	USB ControlBoxx controller	USB, RS-232, electronic interface



## Mechanical drawings



## Plug and Play version provide with:

- RemoteBoxx controller
- Manual shutter
- DC power supply

## Available options

- ACX-HTSK-1 Heat sink for operation in ambient air up to 38°C
- OPT-SHTE-LCX Electro-mechanical shutter
- L1C-MPA for power adjustment
- L1C-AOM for high speed modulation

## Fiber coupling options \*\*

	ACX-FCSMUV or ACX-FCPM	ACX-FCMM-R50
Power stability over 8 hours and within $\pm 3k$	$\pm 2.0\%$	$\pm 2.0\%$
Fiber optical characteristics	3.5 $\mu\text{m}$ core diameter 0.11 N.A	50 $\mu\text{m}$ core diameter 0.22 N.A
Polarization extinction ratio (ACX-FCPM only)	$> 1000 : 1$	n/a
Fiber output connector	FC/APC	FC/APC
Fiber length	2 m	2 m

Warranty : 18 months from shipment date

\*Specifications at nominal power

\*\* Variant are available on demand



## CONTACT US

Oxxius S.A  
4 rue Louis de Broglie  
F-22300 Lannion, France

Phone : +33 296 48 70 28  
sales@oxxius.com  
F-22300 Lannion, France www.oxxius.com